

# Tracor Analytic

Tracor Analytic, Inc.  
1842 Brunnel Drive  
Elk Grove Village, Illinois 60007  
(312) 364 9100

We are pleased to report that the last time we monitored your A/C line (February 28, 1983 through March 7, 1983) we found it to be fully within our factory specifications. The last five days of this test, using the latest Sine Control device, showed no impulse signals (spikes) and no line surges. During that period there were 288 instances of the line sagging. However, the lowest voltage recorded was only 106 VAC, occurring for an average 7.5 cycles. As our specification calls for 115 VAC  $\pm$  10 VAC this is within our recommended operating range.

Dr. Farber  
Clinical Pathology Department, Room 4705  
Cedar Sinai Medical Center  
8700 Beverly Blvd.  
Los Angeles, CA 90048